



Bathing Care

Easy access bathing from
Enjoy Life Mobility.
because we care about
your independence



ELM Enjoy Life Mobility

Enjoy Life Mobility is a well established mobility company and represents the interests of Wibtrac in Spain. Wibtrac is one of the leading manufacturers of Walk in Baths and along with Elm has many years of experience successfully fulfilling the requirements of clients in Europe, Australia, the USA and Canada.

*We are now proud to offer you the Mayfair & Regency range of **Walk in Baths** at exceptional value for money.*

*All our **walk in baths** are manufactured in fibreglass, using the latest technology and materials offering UK Design and UK Quality. All of our **walk in baths** are held at our modern warehouse offering the latest logistical distribution services which allows us to give you a competitive delivery cost and service to wherever you want the bath delivered.*

Our service is based on being honest, reliable with an unprecedented knowledge of our products.

We believe it is imperative to identify all of our customers' needs and objectives and then provide them with an expedient service, giving them professional advice with as much downloadable information as possible.

*We believe that if we are careful and thorough with each and every one of our customers we will ensure that their bathing experience will return as will their quality of life. When you purchase one of our range of Wibtrac **Walk in Baths**.*

***We truly believe** you are buying the best Walk in Bath in the world. All of our baths are hand-made. Every Wibtrac Walk in Bath is Five-Star Factory Tested to ensure that your **Walk in Bath** is in perfect condition when it arrives.*

 *Enjoy Life Mobility*

Bathing safety can be a concern for older adults. Thirty percent of people over the age of 50 experience difficulty getting into or out of their bath or shower, even when safety equipment is installed. Bathing is an essential part of life it cleanses helps lower stress and has a beneficial effect on many ailments. Elm's range of Walk in Baths allow you to once again enjoy a relaxing soak while knowing your safety has been

The decision was made for me when I slipped in the shower, I was so lucky I escaped a serious injury. I now enjoy a regular relax in my new Regency walk in bath Thanks to Enjoy Life Mobility. Winifred Allen



- FACT: You won't beat our prices.
- FACT: You won't get a better made bath.
- FACT: The only bath on the market with stainless steel frame.
- FACT: A lifetime guarantee on the door seal.
- FACT: A safer bathing experience.
- FACT: You are dealing with people who have traded in Spain for 10+ years.



0044 77 666 90 676 or email: baths@enjoylifemobility.com
We guarantee our prices are unbeatable



ELM Enjoy Life Mobility

The Regency

Bathing Care

The Regency range has a Stainless Steel Frame, a foot operated Pop up Waste so you do not have to bend down to insert the plug. The Regency comes in a classic satin white gloss finish. Simply sit down, close the door and when you are safely seated fill the bath to your desired level. Fill it to just above the ankles for that soothing foot massage, above the knees to ease the joints or to the top and ENJOY.

- Bath Height 1041mm - 41 inches
- Bath Width 660mm - 26 inches
- Bath Length 1310mm - 51.5 inches
- Step Height 150mm - 7 inches
- Seat Height 445mm - 17.5 inches
- Seat width 470mm - 18.5 inches



Ask about our Hydro/Spa options



ELM Enjoy Life Mobility

The Mayfair

At last, an affordable Walk-in-Bath all the family can use. It's the same size as conventional full-length baths but has the advantage of the Walk-in-Facility provided by a door at one end. At 1700mm or 1500mm long and with a door either side it makes the Mayfair adaptable to most bathrooms.

- Length 1700mm (66.5") or 1500mm (59.5")
- Width 699mm (27.5")
- Height 571mm (22.5")
- Entry Step Height 175mm (7")
- Door Inward Opening



Easy access bathing from
Enjoy Life Mobility. because
we care about your
independence



Ask about our
Hydro/Spa
options

ELM Enjoy Life Mobility

Thank you for downloading our brochure. We hope you have enjoyed your visit to our website and found it informative. We are always looking to improve our product range and are in the process of introducing a fantastic shower cubical. Please re-visit our site occasionally to see our new products. Produced by Enjoy Life Mobility



ELM Enjoy Life Mobility

Enjoy Life! - Enjoy Life Mobility - Bathing Care

To order your bath
please phone: 0044 776 669 0676
or email us at
baths@enjoylifemobility.com

If you have any questions, have a look at our frequently asked questions on our website at
www.enjoylifemobility.com

or just send us an email and we will answer any queries you may have.

Thank you for you time

ELM Enjoy Life Mobility